

VOICE CHEAT SHEET

To Start: Click the AI button



Denti AI

POCKET DEPTH: In Triplets (434)(323)(212)

NAVIGATION COMMANDS

- **JUMP:** (Command + Tooth) Jump 15/Jump 25 Facial/Jump 31 Lingual
- **UNDO:** Remove last complete command
- **REPEAT:** Repeat multiples of triplets on teeth, (323) Repeat 8 (Up to Repeat 32)
- **PAUSE:** Pause Voice in place until it is Resumed
- **RESUME:** Resume Voice in exact location after it has been Paused

GINGIVAL MARGIN: In Triplets (010)(111)(020)

FIRST SAY: Jump Gingival Margin or Jump Recession to move the cursor to the GM Row

USE SAME COMMANDS FOR PD & GM!

NOTE: PD & GM ARE ABSOLUTELY DEPENDENT ON WHERE THE CURSOR IS LOCATED. YOU MUST BE IN THE PD OR GM ROW TO CHART CORRECTLY

BLEEDING & SUPPURATION: Command/Tooth/Location

- BLEEDING 13 Interproximal
- BLEEDING Quadrant (1-4)
- BLEEDING Sextant (1-6)
- BLEEDING ALL

- SUPPURATION 13 Interproximal
- SUPPURATION Quadrant (1-4)
- SUPPURATION Sextant (1-6)
- SUPPURATION ALL

MOBILITY: Mobility/Tooth/Class or Grade

- MOBILITY 11 Class 1
- MOBILITY 18 Grade 2

FURCATION: Furcation/Tooth + Surface/Class or Grade

- FURCATION 16 ML Class 2
- FURCATION 31 B Grade 1

PLAQUE: Plaque/Location/Light, Medium, Severe

- PLAQUE 31 Severe
- PLAQUE ALL Light

NOTE: BLEEDING/SUPPURATION/MOBILITY/FURCATION/ PLAQUE DO NOT MOVE THE CURSOR. WITHIN THE COMMAND YOU TELL THE AI THE PROPER PLACEMENT

****Finish with "SAVE CHART" Command****

Contact Customer Success:

- support@dentia.ai
- (865) 551-9304